



▶▶ Since 1992

Be healthy for life

Your reliable supplier in Ob&Gyn products ♥



- Fetal/Maternal Monitor
- Fetal Monitor
- Fetal Heart Rate Doppler
- Central Monitoring System






Be Healthy For Life

Sunray SRF618 B6 Fetal Monitor

B6 is a new design with foldable color screen and one shuttle button, which enhances users functionality. Its professional CTG analysis function meets the needs of obstetrical department in offices, clinics and hospitals.

B6 can monitor two pregnant women simultaneously with two-bed CTG graph display.

**Two
in
One**



1200 hours CTG graph storage.

Insight Software / USB



Insight software and USB for data transmission to PC.



Touch key and Shuttle Button



The touch key and shuttle button can help the operator simplify the workflow and work with one button only.



Integrated Wireless Waterproof Probes ▶▶

Integrated wireless waterproof probe module. The probes can be used universally in all of Sunray fetal monitors.

Highly sensitive 12-crystal probes provide wider beam and more stable and homogeneous signal.



Wireless transducer with patent technology
Wide beam and high sensitive waterproof transducer



Upgraded Integrated design

Signals Overlap Verification To Separate Twins And Triplets FHR ▶▶



■ Advanced FHR detection technology

CTG Analysis ▶▶

Professional CTG analysis function meets the needs of obstetrical department in offices, clinics and hospitals.

Real-time analysis function: automatically calculate and analyze the parameter indexes.

Analysis result			
NST Result	Single	Twins	Triplets
Starting Time	0'0"		
Analysis Time	60'0"		
Baseline	144	130	162
Amplitude	15	15	16
Acc Time	22	13	0
Acc Value	22	20	0
FM Times	18	18	18
Acc Times	11	3	0
ED Times	0	0	0
LD Times	0	0	1
VD Times	0	0	0
DEC Times	0	0	1
STVoss	6.3	8.4	4.7

press dial or "menu" to exit

Analysis done Exit

Wired & Wireless Network ▶▶

Wired & wireless network module to connect central monitoring station.

Foldable Color Screen ▶▶

10.2" TFT foldable color screen. Black, green and pink color display interface can be easily switched.

Standard Configuration: FHR, TOCO, FM, AFM

Optional Configuration: Twins FHR, Triplet FHR, Wireless transducer, Touch screen, Li-ion Battery






●● Be Healthy For Life

Sunray **SRF618 K9** Fetal & Maternal Monitor

Monitor The Maternal And Fetal Parameters Simultaneously >>

Monitor the maternal and fetal parameters simultaneously, Fetal parameters: FHR, TOCO and FM; Maternal parameters: NIBP, SPO2, PR, ECG, RESP and TEMP
Both fetal and maternal parameters can be output and printed.

Signals Overlap Verification To Separate Twins And Triplets FHR >>



800 hours CTG graph and maternal parameters storage.

Wired & Wireless Network >>

Wired & wireless network module to connect central monitoring station.

Foldable Color Screen >>

12.1" TFT-LCD foldable color screen. Black, green and pink color display interface can be easily switched.

Standard Configuration: FHR, TOCO, FM, AFM, NIBP, SPO2, PR, ECG, RESP and TEMP

Optional Configuration: Twins FHR, Triplet FHR, Wireless transducer, Touch screen, Li-ion Battery

Integrated Wireless Waterproof Probes ▶▶

Integrated wireless waterproof transducer module. The probes can be used universally in all of Sunray fetal monitors.

Highly sensitive 12-crystal probes provide wider beam and more stable and homogeneous signal.



Wireless transducer with patent technology



Upgraded Integrated design

Analysis result			
NST Result	Single	Twins	Triplets
Starting Time	0'0"		
Analysis Time	00'0"		
Baseline	144	130	162
Amplitude	15	15	16
Acc Time	22	13	0
Acc Value	22	20	0
FM Times	18	18	18
Acc Times	11	3	0
ED Times	0	0	0
LD Times	0	0	1
VD Times	0	0	0
DEC Times	0	0	1
STV(ms)	6.3	8.4	4.7

press dial or "menu" to exit

Analysis done Exit

CTG Analysis ▶▶

Professional CTG analysis function meets the needs of obstetrical department in offices, clinics and hospitals. Real-time analysis function: automatically calculate and analyze the parameter indexes.

Two-bed Monitoring Function ▶▶

K9 can monitor two pregnant women simultaneously with two-bed CTG graph display.

Insight Software / USB ▶▶

Insight software and USB for data transmission to PC.



Touch key and Shuttle Button ▶▶

The touch key and shuttle button can help the operator simplify the workflow and work with one button only.





●●● Be Healthy For Life

Sunray

SRF618 B5 Fetal Monitor

Portable 8.4" TFT color LCD display, 70 degree foldable fetal monitor ▶▶

- ◆ FHR and TOCO display in curves and digitals
- ◆ Waterproof FHR and TOCO probes. Highly sensitive 12-crystal probes provide wider beam, more stable and homogeneous signal
- ◆ Audible and visual alarm of probes off-position, abnormal FHR and unstable signal
- ◆ Built-in wired network connection capability
- ◆ 4 types of self-defined event marks
- ◆ FDA (510K) approved

Standard Configuration: Twins FHR, TOCO, FM, AFM ▶▶

Optional Configuration: wireless network connection ▶▶



Sunray

SRF618 B3 Fetal Monitor

Compact 4.7" smart fetal monitor/flexible workstation ▶▶



SRF618 B3



Workstation



SRF618 B3



SRF618 B3



SRF618 B3

- ◆ Wireless transducers for more comfortable monitoring, FHR and TOCO display in curves and digitals
- ◆ Wireless connection to powerful workstation for data storage, CTG analysis, report printing
- ◆ Audible and visual alarm of probes off-position, abnormal FHR and unstable signal
- ◆ Flexible workstation fits for labor room and obstetric office, up to 9 B3s connection to workstation

Standard Configuration : FHR1, TOCO, FM, AFM ▶▶

Optional Configuration : FHR2 ▶▶





●●● Be Healthy For Life



Sunray SRF618 E Fetal Doppler

Compact type

- ◆ Accurate FHR based detecting based on auto-gain control function
- ◆ Multi-detecting modes of FHR: transient FHR, mean FHR and manual FHR
- ◆ Intellectual cancellation of erroneous FHR
- ◆ Audible alarming on FHR value over-limitation
- ◆ Waterproof USB transducer, easy for replacement
- ◆ Low battery alarming and real-time battery status indicator

Optional Configuration: Earphone ▶▶

Sunray SRF618 D Fetal Doppler

Classic type

- ◆ Accurate FHR detecting with digital back-light LCD display
- ◆ Intellectual cancellation of erroneous FHR and transient FHR displaying
- ◆ Rechargeable battery ensures continuous working more than 14 hours
- ◆ Waterproof transducer with IP68 grade
- ◆ Automatically power off within 2 minutes
- ◆ Low battery alarming and real-time battery status indicator



Optional Configuration: Earphone ▶▶



Be Healthy For Life


Sunray SRF618 S Central Monitoring System



- ◆ Single or multi display windows
- ◆ Wired & wireless central monitoring system solution: Ethernet network, RS485, medical RF, WIFI
- ◆ Patient data analysis function, supporting Krebs, Fischer, NST and CST analysis
- ◆ CTG cases can be preserved, reviewed, commented and printed
- ◆ Connection of up to 200 bedside fetal/maternal monitors, can be extended to 600 monitors

The best total solution for fetal/maternal monitoring

- ◆ Web viewer: supporting doctors with remote monitoring from their own PC, improving communication efficiency
- ◆ Mobile monitoring: view the monitoring information through Pad and cell phone
- ◆ Connect to Hospital Information System(HIS) perfectly
- ◆ Three layers of structure design: database, server, client
- ◆ Basic and advanced alerting immediately informs unforeseen critical events
- ◆ Support statistical evaluation and reporting for management and research
- ◆ Analysis of acceleration, deceleration, long-term variability and short-term variability of FHR
- ◆ Can be connected to Oxford/GE/Philips models





▶ Since 1992

Sunray Medical Apparatus Co.,Ltd.

Manufacturer center:

No.38, Gaoke Rd Gaotang Industry District, Guangshan Er Rd, Guangzhou, China

Marketing center:

4/F, No. 242, Tianhe Dong Road, Guangzhou, P.R. China, 510620

Tel: +86-20-8759-7231 E-mail: overseas@sunray.cn

All rights reserved. Features and specifications are subject to change without notice.

www.sunray.cn