



nu med technologies

the utmost degree of confidence in quality & service

General Surgery

Discover “numed” make – A absolute Prerequisite key for quality of its range of Surgical Instruments



“numed” Instruments for General Surgery, surgical Oncology, Micro Vascular, Plastic Surgery, Urology, Gynaecology

A-202, II FLOOR, CHARAN SINGH COMPLEX,
M.B. ROAD, SAKET, NEW DELHI – 30.
PH. 011- 40506330, FAX. 40506331
Mob.- 09818274945, 09910222245
EMAIL : numedtech@gmail.com



Total Quality Assurance

As a medical professional, you want to ensure the utmost quality and integrity of your surgical instruments at all times. You also want to protect your considerable investment in these high-grade medical devices.

USP OF “nu med” Surgical Instruments

Steel Components:

We are currently using AISI 410, AISI 420, and 18/8 type of stainless steel and according to the items specifications.

Anti Rust Protection

We have accurate furnaces to obtain hardness of C-40 to C-45 on Forceps and between C-50 TO C-55 on Scissors and Sharpened Instruments.

Instruments Quality Warranty:

We offers 2 year's quality warranty against rust, manufacturing defect. Free replacement is offered in case of any manufacturing defect, though there is a minor chance of such happening.

Specialty Products:

“nu med” special products are carefully designed according to end user application to enrich surgical techniques.



No. 3 Handle
BP-000-03



No. 4 Handle
BP-000-04



No. 7 Handle
BP-000-07

DT-002-35 HUMBY
DERMATOME FOR
UNIFORM GRAFTING



1/2



No. 3L Handle
BP-000-03L




No. 3L Handle
BP-000-03A



No. 3L Handle
BP-000-04A



DT-001-19 Silver Razor
Dermatome



Backhaus Towel Clamp
 CMP-002-12 – 4" 12 cm
 CMP-002-13 – 5" 13 cm
 CMP-002-15 – 6" 15 cm



Jones Martin Towel Clamp
 CMP-001-09 – 2 3/4" 9 cm
 CMP-001-11 – 3 1/2" 11 cm



TFC-003-11 Schadel
 Towel Clamp



TFC-004-11 Jones Martin
 Towel Clamp 8.5 cm



TFC-005-15 Roder
 Towel Clamp
 6", 15 cm



TFC-006-16
 Bernhard Towel
 Clamp 16.5 cm



CMP-003-20 Lane
 Towel Clamp 8" 20 cm



Tubing Clamps with Safety Guard
 CMP-004-15 – 6", 15 cm
 CMP-004-18 – 7", 18 cm
 CMP-004-20 – 8", 20 cm



Tubing Clamps with Safety Guard
 CMP-004-15 – 6", 15 cm
 CMP-004-18 – 7", 18 cm
 CMP-004-20 – 8", 20 cm




RO-VS-18 GUNNAR-HEY
 Tube Milking Roller 7 1/4",
 18.5 cm




Lister Sinus Forceps
 FR-001-13 – 5" 13 cm
 FR-001-15 – 6" 15 cm



ULRICH Sponge & Dressing Forceps
 FR-002A-20 – 8", 20 cm
 FR-002A-25 – 10", 25 cm




RAMPLEY Sponge Holding Forceps
 FR-002B-18 – 7", 18 cm
 FR-002B-25 – 10" 25 cm



Foerster Sponge Holding Forceps
 FR-002-18 – 7", 18 cm
 FR-002-20 – 8", 20 cm
 FR-002-25 – 10", 25 cm





AS-002-20 Poole
Suction 10 MM ST



AS-002-20 Poole
Suction 10 MM Cvd



AS-005-06 – 6 Fr.
AS-005-08 – 8 Fr.
AS-005-10 – 10 Fr.
AS-005-12 – 12 Fr.




Frazier Suction St & Cvd Pattern



AS-001-20 – 8” 20 cm
AS-001-25 – 10” 25 cm
Magill Forceps



YANKAUER SUCTION
CANULLA 11”



Director Probe & Buttoned Probe

“numed” SURGICAL SCISSORS



Iris Delicate Scissor ST
 SR-024-11 - 4 1/2" 11.5 cm
 SR-024-14 - 5 1/4" 14 cm





Iris Delicate Scissor CVD
 SR-025-11 - 4 1/2" 11.5 cm
 SR-025-14 - 5 1/4" 14 cm



Baby Metzenbaum ST
 SR-010-11 - 4 1/2" 11.5 cm



Baby Metzenbaum CVD
 SR-010-11 - 4 1/2" 11.5 cm



Standard Surgical Scissors SH/BL

STRIGHT PATTERN SH/BL

- SR-001-15 - 6" 15 cm
- SR-001-18 - 7" 18 cm
- SR-001-20 - 8" 20 cm

CURVED PATTERN SH/BL

- SR-002-15 - 6" 15 cm
- SR-002-18 - 7" 18 cm
- SR-002-20 - 8" 20 cm


STRIGHT PATTERN SH/SH

- SR-003-15 - 6" 15 cm
- SR-003-18 - 7" 18 cm
- SR-003-20 - 8" 20 cm



**MAYO DISSECTING SCISSORS
STRIGHT**

- SR-004-13 – 5", 13, cm
- SR-004-15 – 6", 15 cm
- SR-004-18 – 7", 18 cm
- SR-004-20 – 8", 20 cm
- SR-004-23 – 9", 23 cm




**MAYO DISSECTING SCISSORS
CURVED**

- SR-005-13 – 5", 13, cm
- SR-005-15 – 6", 15 cm
- SR-005-18 – 7", 18 cm
- SR-005-20 – 8", 20 cm
- SR-005-23 – 9", 23 cm



**METZENBAUM OPERATING
SCISSORS STRIGHT**

- SR-012-15 – 6", 15 cm
- SR-012-18 – 7", 18 cm
- SR-012-20 – 8", 20 cm
- SR-012-23 – 9" 23 cm
- SR-012-25 – 10" 25 cm



**METZENBAUM OPERATING
SCISSORS CURVED**

- SR-013-15 – 6", 15 cm
- SR-013-18 – 7", 18 cm
- SR-013-20 – 8", 20 cm
- SR-013-23 – 9" 23 cm
- SR-013-25 – 10" 25 cm



**METZENBAUM OPERATING
SCISSORS STANDARD**



METZENBAUM FINO METZENBAUM FINO
 SR-014-15 – 6” SR-014-15 – 6”
 SR-014-18 – 7” SR-014-18 – 7”
 SR-014-20 – 8” SR-014-20 – 8”

KELLY (FISTULA) SCISSORS
 SR-026-16 – ST
 SR-027-16 – CVD

Steven Delicate Surgical Scissors Fine
 SR-028-11 – St
 SR-029-11 – Cvd




Lister Bandage Scissor
 6” – SR-047-15
 8” – SR-047-20

LORENZ Bandage Scissor 23 cm

BURNS Plaster Scissor 24 cm

CASTROVIEJO Micro-Spring type/ Vascular Scissors Flat Handle
A - 7" 18 cm & B - 8" 20 cm (Spring Pattern Available in TC Also.)



25 DEGREE ANG
SR-052-25

45 DEGREE ANG
SR-052-45

60 DEGREE ANG
SR-052-60

90 DEGREE ANG
SR-052-90


125 DEGREE ANG
SR-052-125



5" St – SR-049-13
 6" St – SR-049-15
 7" St – SR-049-18
 8" St – SR-049-20



5" Cvd – SR-048-13
 6" Cvd – SR-048-15
 7" Cvd – SR-048-18
 8" Cvd – SR-048-20



Potts Smith Vascular Scissors 7 1/2" 19 cm
Ring Handle
 SR-038-25 – 25 Deg
 SR-038-45 – 45 Deg
 SR-038-60 – 60 Deg
 SR-040-125 – 125 Deg REVERSE CUTTING



Littauer Ligatures Scissors
5" – SR-037-13




Heath Ligatures Scissors
6 1/4" – SR-036-15




Spencer Scissors
5 1/2" – SR-035-14




Epistomy Scissors
5 1/2" – SR-045-15



SR-044-20 SIMS Uterine St
Gynecological Scissors



SR-044-20 SIMS Uterine Cvd
Gynecological Scissors



UNIVERSAL Wire Cutter
For soft wire 1.2 mm

1 Blade Serrated

“numed” Tungsten Carbide Scissors



Isis Scissor Serrated
4 1/2" St - TSR-024-11
5 1/4" St - TSR-024-14



Isis Scissor Serrated
4 1/2" St - TSR-025-11
5 1/4" St - TSR-025-14




Iris Scissor 4 1/2" Angle
on side TSR-025-11 A



Epistomy Scissors
SR-045-15 5 1/2"



DEAN Gum scissors
SR-033-17 Cvd 6 3/4"



TC KELLY Scissors
St - TSR-026-16
Cvd - TSR-027-16



TC Kilner Fine Scissors
St - TSR-016-12
Cvd - TSR-017-12



1 Blade Serrated



TC UNIVERSAL Wire Cutter
For soft wire 1.2 mm

Mayo Stille Dissecting Scissors with Tungsten Insert



Stright Pattern

6" St - TSR-006-15

7" St - TSR-006-18

Curved Pattern

6" Cvd - TSR-007-15

7" Cvd - TSR-007-18

Mayo Dissecting Scissors with Tungsten Insert



Stright Pattern

6" St - TSR-004-15

7" St - TSR-004-18

9" St - TSR-004-23


Curved Pattern

6" Cvd - TSR-005-15

7" Cvd - TSR-005-18

9" Cvd - TSR-005-18

**Metzenbaum Operating Scissors with Tungsten Insert
Standard Pattern**



Stright Pattern

6" St - TSR-012-15

7" St - TSR-012-18

8" St - TSR-012-20

9" St - TSR-012-23

10" St - TSR-012-25

Curved Pattern

6" Cvd - TSR-013-15

7" Cvd - TSR-013-18

8" Cvd - TSR-013-20

9" Cvd - TSR-013-23

10" Cvd - TSR-013-25



TSR-039-25 Potts Scissor 25 Deg



TSR-039-45 Potts Scissor 45 Deg



TSR-039-60 Potts Scissor 60 Deg



TSR-039-90 Potts Scissor 60 Deg



TSR-039-125 Reverse Potts Scissor 125 Deg



“numed” Traditional Micro Scissors

Ideal for cutting thick vascular and scar tissue
 Round or flat handle
 Lengths of 20, 22, 24, 26 cm
 Straight or curved tips




“numed” Ultra Fine Micro Scissors

Extremely thin blades , gold springs
 Ideal for dissection of arachnoid
 Round or flat handle
 Lengths of 20, 22, 24, 26 cm
 Straight or curved tips



“nuned” Meningioma/Tentorium Micro Scissors
Designed for opening the tentorium
Short blades with sharp tips, black springs
Lengths of 20, 22 cm Flat handle
Straight or curved tips



Titanium Micro Scissor
Round Handled
TSR-053-08 – 3 ¼” St
TSR-054-08 – 3 ¼” Cvd



Castroviejo Style Micro Scissors
Round or Flat Handle
18cm Length
Curved Pattern
TSR-048-15 – 6” Cvd Blade
TSR-048-18 – 7” Cvd Blade
TSR-048-20 – 8” Cvd Blade
Straight Pattern
TSR-049-15 – 6” St Blade
TSR-049-28 – 7” St Blade
TSR-049-20 – 8” St Blade



Barraquer
Micro Needle Holder St Barraquer
Micro Needle Holder Cvd

Straight Pattern:

MNH-001-11 - 4 1/2" 11 cm
MNH-001-14 - 5" 14 cm
MNH-001-15 - 6" 15 cm

Curved Pattern:

MNH-002-11 - 4 1/2" 11 cm
MNH-002-14 - 5" 14 cm
MNH-002-15 - 6" 15 cm



Castroviejo Micro Needle Holder
Short Handle 5 1/4" 13 cm



Castroviejo Micro Needle Holder
with Lock Flat Handle Grip

Straight Pattern:


MNH-004-14 - 5" 14 cm
MNH-004-18 - 7" 18 cm
MNH-004-20 - 8" 20 cm
MNH-004-25 - 10" 25 cm

Curved Pattern:

MNH-005-18 - 7" 18 cm
MNH-005-20 - 8" 20 cm



Kalt Needle Holders
5 1/2" 14 cm



HALSEY Needle Holder
13 cm 5"




BAUMGARTNER
Needle Holder
13 cm 5"




Crile Wood Needle Holder

NH-003-15 – 6” 15 cm
 NH-003-18 – 7” 18 cm
 NH-003-20 – 20” 20 cm




NEWORLEANS Needle Holder
 7”, 18 cm



Mayo Needle Holder
 Standard Pattern

NH-001-15 – 6” 15 cm
 NH-001-18 – 7” 18 cm
 NH-001-20 – 8” 20 cm
 NH-001-25 – 10” 25 cm



WANGENSTEEN Needle Holder
 27 cm 10 1/2”



Mayo Needle Holder
 Delicate Pattern

NH-002-15 – 6” 15 cm
 NH-002-18 – 7” 18 cm
 NH-002-20 – 8” 20 cm
 NH-002-25 – 10” 25 cm

RYDER NEEDLE HOLDER

NH-004-15
 NH-004-18
 NH-004-20






NH-005-21 Heany
Needle Holder 21 cm 8 1/4"

NH-008-18 Olsen-Hegar
Needle Holder 7"




NH-006-25 Finochietto
Needle Holder 25 cm 10"

Gillis Needle Holder
15 cm, 6"



NH-007-23 Stratte Needle Holder
24cm 9", Curved on side

Jhonson Needle Holder
27cm, 10 1/2"





TC MAYO-Hegar
STANDARD Pattern

TC MAYO-Hegar
STANDARD


- TNH-001-15
- TNH-001-18
- TNH-001-20
- TNH-001-15



TC MAYO-Hegar-Delicate/
Hegar Light Pattern

TC MAYO-Hegar-Delicate/
Hegar Vascular


- TNH-002-15
- TNH-002-18
- TNH-002-20
- TNH-002-15



TC DE BAKEY Fine

TC DE BAKEY Fine Needle Holder

- TNH-005-15 – 6" 15 cm
- TNH-005-18 – 7" 18 cm
- TNH-005-20 – 8" 20 cm



TC RYDER VASCULAR

TC RYDER VASCULAR
Needle Holder Dotted Grip
& Smooth Grip

- TNH-004-15 – 6" 15 cm
- TNH-004-18 – 7" 18 cm
- TNH-004-20 – 8" 20 cm




TC GILLIES
Needle Holder 6"



TC NEWORLEANS
Needle Holder 7"



TC SAROT
Needle Holder 7" & 8"



TC CARROL
Needle Holder 6 1/4" 16 cm



TC CRILE-WOOD Needle Holder
TNH-003-15 – 6" 15 cm
TNH-003-18 – 7" 18 cm
TNH-003-20 – 8" 20 cm



TC WANGENSTEEN
10 1/2" 26 cm



TC HEGAR VASCULAR
Needle Holder 8" – 20 cm
& 9 1/2" – 24 cm



TC Dietrich
Needle Holder, 7" & 8"



TC Parrel
Needle Holder, 7"



TC Hegar-Baumgartner
TNH-001-13



TC Mathieu
Needle Holder, 6" & 8"



TC Heany
Needle Holder, 8"



TC Bozemann
Needle Holder, 8" & 10"



TC Olsen-Hegar
Needle Holder, 7" & 8"



TC Lichtenberg Ryder
Needle Holder, 6 1/2" & 8"




BARRAQUER Micro
Needle Holder
MNH-001-11 - 5½"



BARRAQUER Micro
Needle Holder
MNH-001-11 - 5½"



TC KALT 14 cm 5 ½"
Needle Holder



Halsey 13 cm 5"
Needle Holder



TMNH-004-18 – 7" 18 cm
TMNH-004-20 – 8" 20 cm
TMNH-004-25 – 10" 25 cm

Castroviejo TC Micro Needle Holder with Lock
Flat Handle Grip



GIANNINI Long Micro
Needle Holder 9 ½"



Castroviejo TC Micro
Needle Holder 13 cm 5"



Micro Jeweler's Forceps
MFP-001-13, 4 3/4" 13 cm



Mini-Adson Dissecting Forceps
FP-006-14 4 3/4" 14 cm
FP-006-15 6" 15 cm

Non Tooth



Micro-Adson Dissecting Forceps
FP-006-14 4 3/4" 14 cm
FP-006-15 6" 15 cm

Non Tooth



Yasargil Micro Tumor and Tissue Grasping Forceps, baynet-Shaped

0,4 mm



Adson Dissecting Forceps
FP-005-14 - 4 3/4" Plain 14 cm
FP-005-14 - 4 3/4" Plain 14 cm

Non Tooth



Yasargil Ring Micro Grasping Forceps baynet-Shaped 7", 8", 10"

3 mm




FP-017-17 GERALD delicate dissecting forceps Plain 6 3/4"

Non Tooth

TAYLOR Delicate Dissecting Forceps
17,5 cm, 7"

FP-010-18 - St
FP-011-18 - Cvd



Non Tooth

1/2


1/2



Non Tooth



Non Tooth



Non Tooth

Potts-Smith delicate forceps

- FP-014-18 7"
- FP-014-20 8"
- FP-014-25 10"
- FP-014-30 11"

McIndoe delicate forceps
FP-009-15, 15 cm 6"

Taylor Cushing Forceps
FP-012-18, 17.5 cm 7"



Non Tooth

WAUGH Delicate Forceps
FP-016-20 8" 20 cm



1/4

Brophy Delicate Forceps
8" 20 cm St & Cvd



1/4



Non Tooth

GILLIES
Delicate Forceps
6" 15 cm



Non Tooth

MOD-USA Dissecting Forceps
Standard Pattern

- FP-019-14 5" 14 cm
- FP-019-15 5 3/4" 14.5 cm
- FP-019-16 6 1/4" 16 cm
- FP-019-18 7" 18 cm
- FP-019-20 8" 20 cm
- FP-019-25 10" 25 cm


MOD-USA Dissecting Forceps
Narrow Pattern

- FP-021-14 5" 14 cm
- FP-021-15 5 3/4" 14.5 cm
- FP-021-16 6 1/4" 16 cm
- FP-021-18 7" 18 cm
- FP-021-20 8" 20 cm
- FP-021-25 10" 25 cm



Non Tooth

1/1




Adson Vascular



Adson Forceps
 FP-007-14 4 3/4" 14 cm
 FP-007-15 6" 15 cm



Micro-Adson Forceps
 FP-008-14 4 3/4" 14 cm
 FP-008-15 6" 15 cm




Micro-Adson Forceps
 FP-008-14 4 3/4" 14 cm
 FP-008-15 6" 15 cm




Adson Brown
 Forceps 9 x 9 Teeth



Cushing Delicate
 Forceps 6 3/4" 17 cm



Taylor-Cushing
 Forceps 1 x 2 Teeth
 17.5 cm, 7"



Potts-Smith
 FP-015-18 7"
 FP-015-20 8"
 FP-015-25 10"
 FP-015-30 11"



WAUGH
 Delicate Forceps
 FP-016-20



GERALD
 Delicate Forceps
 FP-018-17



GILLIS
 Delicate Forceps
 FP-018-17



Tooth

MOD-USA Dissecting Forceps
Standard Pattern, 1 x 2 Tooth

FP-020-14	5"	14 cm
FP-020-15	5 3/4"	14.5 cm
FP-020-16	6 1/4"	16 cm
FP-020-18	7"	18 cm
FP-020-20	8"	20 cm
FP-020-25	10"	25 cm

MOD-USA Dissecting Forceps
Narrow Pattern, 1 x 2 Tooth

FP-022-14	5"	14 cm
FP-022-15	5 3/4"	14.5 cm
FP-022-16	6 1/4"	16 cm
FP-022-18	7"	18 cm
FP-022-20	8"	20 cm
FP-022-25	10"	25 cm



Tooth



Tissue Forceps Cvd
5 1/2" 14 cm



LANE Forceps 3 x 2
5 1/2" 14 cm
8" 20 cm




BONNEY Forceps
Broad Tip 7" 18 cm



ADLERKREUTZ Vascular
& Tumor Forceps 5" 14 cm
& 8" 20 cm




SINGLEY-TUTTLE Lung
Forceps 7" 18 cm



RUSSIAN Tissue
Forceps 20 , cm 8 "



COLLIN-DUVAL Intestinal
Forceps 20 , cm 8 "



TC Micro Adson Dissecting

FP-005-14 – 4 3/4" 14 cm
 FP-005-15 – 6" 15 cm




TC Adson Dissecting

FP-005-14 – 4 3/4" 14 cm
 FP-005-15 – 6" 15 cm



TC MINI Adson Dissecting

FP-005-14 – 4 3/4" 14 cm
 FP-005-15 – 6" 15 cm



TC DeBakey Vascular Dissecting Forceps with Autramatic Serration

FP-026-15 – 6"
 FP-026-20 – 8"
 FP-026-25 – 10"



Tooth



Tooth



Tooth

TC Micro Adson Dissecting
TFP-007-14 – 4 3/4" 14 cm
TFP-007-15 – 6" 15 cm

TC Micro Adson Dissecting
TFP-007-14 – 4 3/4" 14 cm
TFP-007-15 – 6" 15 cm

TC ADSON BROWN
Forceps 9 x 9 Teeth



TC Gillis 6"
Plain



TC Gillis 6"
Tooth



TC Mc Indoe 6"
15 cm Tip 0.4
mm



Potts Smith
6", 15 cm
7", 18 cm
8", 20 cm



Cushing 7", 18 cm



Gerald 6 3/4"
17.5



TC EVANS 25 cm,
10" 0,4 mm



HALSTED Micro-Mosquito
Forceps ST
FR-003-12 – 4" 12 cm
FR-003-13 – 5" 13 cm



HALSTED Micro-Mosquito
Forceps CVD
FR-004-12 – 4" 12 cm
FR-004-13 – 5" 13 cm



Dandy Forceps
FR-009-14 – 5 1/2" 14cm




Kelly Haemostatic
Forceps ST
FR-007-15 – 6" 15cm
FR-007-18 – 7" 18cm
FR-007-20 – 8" 20cm


Kelly Haemostatic
Forceps CVD
FR-008-15 – 6" 15cm
FR-008-18 – 7" 18cm
FR-008-20 – 8" 20cm



Crile Haemostatic Forceps
FR-010-14 - 5 1/2" 14cm
FR-010-15 – 6" 15cm

Crile Haemostatic Forceps
FR-010-14 - 5 1/2" 14cm
FR-010-15 – 6" 15cm






Roberts Ligature Forceps
FR-019-23 9" ST 24 cm
FR-020-23 9" CVD 24 cm



FR-022-23
Moynihan Forceps
Thick Tips



Rochester-Pean ST



Rochester-Pean ST
FR-012-15 – 6" 15 cm
FR-012-18 – 7" 18 cm
FR-012-20 – 8" 20 cm
FR-012-25 – 10" 25 cm

Rochester-Pean CVD
FR-013-15 – 6" 15 cm
FR-013-18 – 7" 18 cm
FR-013-20 – 8" 20 cm
FR-013-25 – 10" 25 cm



Rochester-Pean CVD

Spencerwell Forceps ST
FR-005-15 – 6" 15 cm
FR-005-18 – 7" 15 cm
FR-005-20 – 8" 15 cm


Spencerwell Forceps CVD
FR-006-15 – 6" 15 cm
FR-006-18 – 7" 15 cm
FR-006-20 – 8" 15 cm



ROCHSTER
CARMALT ST
FR-013S-15
FR-013S-20



ROCHSTER
CARMALT CVD
FR-013C-15
FR-013C-20



Kocher-Delicate Forceps St

FR-014-13 – 5" 13 cm

FR-014-14 – 5 1/2" 14 cm

Kocher-Oschner Forceps ST

FR-014-15 – 6" 15 cm

FR-014-18 – 7" 18 cm

FR-014-20 – 8" 20 cm

FR-014-25 – 10" 25 cm



Kocher-Delicate Forceps Cvd

FR-015-13 – 5" 13 cm

FR-015-14 – 5 1/2" 14 cm

Kocher-Oschner Forceps Cvd

FR-015-15 – 6" 15 cm

FR-015-18 – 7" 18 cm

FR-015-20 – 8" 20 cm

FR-015-25 – 10" 25 cm



Gemini Ligature Forceps

Slightly Curved Tips

FR-014-15 – 6" 15 cm

FR-014-18 – 7" 18 cm

FR-014-20 – 8" 20 cm



FR-016-14 Baby-Mixer
Ligature Forceps 5 1/2"
Semi Cvd



FR-017-14 Baby-Meeker
Ligature Forceps 5 1/2"
Full Cvd



Mixer Right Ang
Forceps

FR-018-15 – 6" 15 cm

FR-018-20 – 7" 18 cm

FR-018-23 – 8" 20 cm


FR-018-25 – 10" 25 cm




LAHEY Vertical Serration
FR-020-19 – 7 1/2" 19 cm
FR-020-23 – 9" 23 cm



OVERHOLT
MIXTER 21 cm, 8 1/4"



GUYON Kidney
Pedicle Clamp
21 cm, 8 1/4"




GUYON Kidney
Pedicle Clamp
23 cm, 9"



CMP-008-21 WERTHEIM
CULLEN Kidney/Hyst. Clamp
8 1/2" 21 cm



CMP-013-23 WERTHEIM
DE BAKEY Hysterectomy/
Parametrium
Clamp 9 1/4" 23.5 cm



Heany Hysterectomy
Forceps Double Tooth
CMP-007-21 – 8 1/4", 21 cm
CMP-007-23 – 9", 23 cm


Wound & Nerve Hook Retractors



Length 15,5 cm, 6"
 Tip A 0,3 mm 25°
 Tip B 0,6 mm 25°
 Tip C 1,0 mm 25°

Length 15,5 cm, 6"
 Tip A 0,3 mm 45°
 Tip B 0,6 mm 45°
 Tip C 1,0 mm 45°

Length 15,5 cm, 6"
 Tip A 0,3 mm 90°
 Tip B 0,6 mm 90°
 Tip C 1,0 mm 90°



SHAMBAUGH Micro Surgery Hooks

LUCAE 14 cm, 5½" Micro Ear Hooklets

CRILE 14,5 cm, 5¾" & 20 cm, 8"

CUSHING 16 cm, 6½" 2 Prongs & 3 Prongs



Blunt



Sharp




CUSHING Nerve Hook 20 cm, 8"

GRAEFE Strabismus hooks 13 cm, 5"

CRILE 18 cm, 7" Sharp

CUSHING Nerve Hook 17 cm, 6 3/4"



CUSHING Nerve Root Retractor 16 cm, 6 1/2"



Gillis Skin retractor, sharp 18, cm, 7"



Gillis Converse Dura & skin retractor, sharp 18, cm, 7"



KILNER Skin Retractor, Sharp 13 cm, 5"





Fig 1: RT 1 J
Fig 2: RT 2 J

CUSHING Vein Retractors 20 cm, 8"

- Sizes Available
- 8 mm
 - 10 mm
 - 12 mm
 - 14 mm
 - 16 mm
 - 18 mm




JOSEPH Hooks for plastic surgery, sharp 16 cm, 6 1/4"



- Sizes Available
- 8 mm
 - 10 mm
 - 12 mm
 - 14 mm
 - 16 mm
 - 18 mm

DESMARRES Vein Retractors 13 cm, 5"






De-Compression Hook

1/2

CUSHING
Vein Retractors/
20 cm, 8"



1/2

Sizes Available:
10 x 6 mm
17 x 5 mm
20 x 6 mm
22 x 8 mm

MINI LANGENBECK
16 cm, 6 1/4"



1/2

8 x 4 mm


15 x 5,5 mm

RAGNELL-DAVIS 14 cm,
5 1/2" 15 x 5,5 mm



1/2

US-ARMY/PARKER
LANGENBECK
20 cm, 8"



1/4

Sharp

1/4

1/2

SENN-MILLER Retractors,
Double-ended 16 cm, 6 1/4"




1/4

Blunt

1/4

1/2



1/4

1/4

1/2

1/2

VOLKMANN Retractors
11,5 cm, 4 1/2 " Sharp &
Blunt

Sharp 2 Prongs
Sharp 3 Prongs
Sharp 4 Prongs

Blunt 2 Prongs
Blunt 3 Prongs
Blunt 4 Prongs



1/2


RT-037-18 Czerney/ Mathieu, Double-
ended 16 cm, 6 1/4" 40 x 15 mm

VOLKMANN Retractors Sharp Pattern Length 8"




1 Prong

RT-038-01



2 Prongs

RT-038-02




3 Prongs

RT-038-03



4 Prongs

RT-038-04



5 Prongs


RT-038-05

VOLKMANN Retractors Blunt Pattern Length 8"




1 Prong

RT-039-01




2 Prongs

RT-039-02




3 Prongs

RT-039-03



4 Prongs

RT-039-04



5 Prongs


RT-039-05

40 x 40 mm




Langenbeck Retractor
Sizes Available:

- RT-040-01 - 30 x 11 mm
- RT-040-02 - 30 x 14 mm
- RT-040-03 - 30 x 16 mm
- RT-040-04 - 40 x 11 mm
- RT-040-05 - 40 x 15 mm
- RT-040-06 - 50 x 11 mm




Pattern 1




Pattern 1

Langenbeck Retractor 21 cm, 8 1/4" Blunt




45 x 50 mm

Flat Blade
60 x 20 mm



Kocher Langenbeck Retractor
Sizes Available:


- RT-042-01 - 22 x 06mm
- RT-042-02 - 25 x 08mm
- RT-042-03 - 35 x 08mm
- RT-042-04 - 35 x 11mm
- RT-042-05 - 35 x 15mm
- RT-042-06 - 40 x 11mm
- RT-042-07 - 55 x 11mm
- RT-042-08 - 70 x 15mm
- RT-042-09 - 80 x 11mm
- RT-042-10 - 80 x 15mm



ISRAEL Retractor 25 cm, 10" Blunt

LANGENBECK 21,5 cm, 8 1/2"


Kocher Langenbeck 22 cm, 8 1/2"



Doyen Retractor Long Blade
 Sizes Available:

- RT-044-01A - 55 x 35 mm
- RT-044-02A - 55 x 45 mm
- RT-044-03A - 85 x 35 mm
- RT-044-04A - 85 x 45 mm


Doyen Abdominal retractor
 23 cm, 9"



Fritsch/ Doyen Retractor
 Handle Vertical to Blade

- RT-044-01B - 30 x 30 mm
- RT-044-02B - 32 x 40 mm
- RT-044-03B - 40 x 50 mm
- RT-044-03B - 40 x 60 mm
- RT-044-04B - 45 x 75 mm


FRITSCH Abdominal Retractor
 23 cm, 9"



Doyen Abdominal retractor
 25 cm, 10" 50 x 85 mm

RICHARDSON/ KELLY
 26 cm, 10 1/4"

- RT-043-01 - 28 x 20 mm
- RT-043-02 - 36 x 28 mm
- RT-043-03 - 44 x 38 mm
- RT-043-04 - 52 x 22 mm




RICHARDSON-EASTMAN SET
 49 x 38 mm, 28 x 20 mm
 26 cm, 10 1/4"




RT-050-01 ST. MARKS Deep Pelvic Retractor 13" 60 mm x 45 mm FIG 1



RT-050-02 ST. MARKS Deep Pelvic Retractor 11 1/2" 35 mm x 45 mm FIG 3



RT-050-02 ST. MARKS Deep Pelvic Retractor 13" 60 mm x 45 mm FIG 3



RT-047-01 - 127 x 40 mm
RT-047-02 - 127 x 62 mm

Harrington Retractor 12 ½", 32 cm



Davidson Retractor 6 ¼", 16 cm