

Medical Suction Unit



Suction units **numed SU 110** (40l/min or 50 l/min) present a new generation of suction units.

This designates them for wider using in surgery, Gynecology, and thorax draining, at as well as for suction from air passages.

There is no service-needed vacuum pump into

| | |
|---------------------------|---|
| Voltage | 220+/-10% V, 50Hz, Single phase AC |
| Max negative pressure | -600 +/-10mmHg (-80kPa) |
| Pumping Rate | 0.5 HP, Heavy Duty Pump, Oil Free |
| Work/ Usage Environment | Temperature 10 deg C to 50 deg C, RH 10% to 90% |
| Bottle Volume | 1000 ml Polycarbonate Bottles |
| Noise | Less than 50 DB |
| Rust Proof Naylor Castors | Mobile with Handle |
| Certificates | ISO & CE Certified |