




nu med technologies

the utmost degree of confidence in quality & service

Gynecology & Obstetrics

Comprise "nu med" specialty vaginal devices to enrich your practice



A-202, II FLOOR, CHARAN SINGH COMPLEX,
M.B. ROAD, SAKET, NEW DELHI – 30.
PH. 011- 40506330, FAX. 40506331
Mob.- 09818274945, 09910222245
EMAIL : numedtech@gmail.com



Total Quality Assurance

As a medical professional, you want to ensure the utmost quality and integrity of your surgical instruments at all times. You also want to protect your considerable investment in these high-grade medical devices.

USP OF "nu med" Surgical Instruments

Steel Components:

We are currently using AISI 410, AISI 420, and 18/8 type of stainless steel and according to the items specifications.

Anti Rust Protection

We have accurate furnaces to obtain hardness of C-40 to C-45 on Forceps and between C-50 TO C-55 on Scissors and Sharpened Instruments.


Instruments Quality Warranty:

We offers 2 year's quality Warranty against rust, manufacturing defect. Free replacement is offered in case of any manufacturing defect, though there is a minor chance of such happening.

Specialty Products:


"nu med" special products are carefully designed according to end user application to enrich surgical techniques.

Bi-Valve Vaginal Speculums



Cusco Broad Blade

GY-001-01 75 × 32 mm
GY-001-02 85 × 35 mm
GY-001-01 100 × 35 mm



Cusco Narrow Blade

GY-001-01 80 × 22 mm
GY-001-02 100 × 25 mm
GY-001-01 105 × 27 mm




**Cusco Virgin
- Menopause**



Graves Vaginal Speculum

GY-002-01 75x20mm Small
GY-002-02 95x35mm Medium
GY-002-03 115x35mm Large




**COLLIN LANDAU, upper blade of
specula shortened 85 X 30 mm**

LATERAL SINGLE BLADE SPATULA



KRISTELLER GY-003-04
Sims Single Ended 'L' Shape
Speculums Medium Size



DOYEN GY-003-04 Single
Ended 'L' Shape Speculums
55 X 35 (S), 85 X 35 (M),
115 X 36 (L)




AUVARD Speculum with
detachable weight 80 x 38
mm



GY-003-01 70x30, 65x35mm Small
GY-003-02 70x35, 70x30mm Medium
GY-003-03 90x40, 85x35mm Large


Sims Vaginal Speculums Double Ended



**MATHIEU VAGINAL SPATULA &
RETRACTORS**

PIC A – 100 x 14 TO 26 MM

PIC B – 95 x 16 TO 26 MM




Bi-Lateral Self Retaining Vaginal Retractor For Colposcopy




Bi-Lateral Vaginal Retractor Non
Conductor Standard Pattern
Specially For LEEP



GY-023-25 **KOGAN** Endocervix
speculum for cervix and uretra
during colposcopy 24 cm, 9,1/2"



CERVICAL PUNCH BIOPSY



SCHUMACHER Biopsy Forceps suitable for colposcopy, 240mm, 9-1/2"



GY-016-21 **Schubert** Cervical Biopsy Forceps 4.5 x 8mm 10 1/4" Ang



GY-016-21 **Vandoren** Cervical Biopsy Forceps 4.5 x 8mm 10 1/4" Ang



16868 GY IUD FORCEP
INTRA UTERINE DEVICE

Instrument for easy finding of the non perforated IUD FORCEP lost IUD and removal of foreign substance. The malleable shaft of the IUD forceps is well adapted to the anatomic conditions.



GY-020-25 Uterine Depressor/
Wall Retractor

Uterine Probes & Endometrical Biopsy Curette



GY-010-30 UTERINE
SOUND/ PROBE



GY-010-30 UTERINE
SOUND/ GRADUATED




Novak GY-011-23
0.4mm with 4 Spikes



Randall GY-011-23
0.4mm 9" 1 Spike




OB-001-27 Wrigley Outlet Forceps 11"



"numed" Murles Vertis Lower Segment
Caesarian Obstetrical Forceps

Hegar Uterine Dilator Female



GY-004-00 Set of 10 Double
Side calibrated

Uterine Curette



GY-013-00
BLAKE Double Ended Sh/BI
Size 0,1,2,3



GOODEL Uterine Dilator 34 cm
13 1/2"



NH-006-25 FINOCCHIETTO "NDVH"
NEEDLE HOLDER 10" 25 CM


Uterine Tenaculum Forceps



GY-006-25
BRAUN 10" 25 cm



GY-006-25 TEALE
VULSULLUM 10" 25 cm



GY-008-25
SCHROEDER
10" 25 cm




GY-008-25 LAHEY
10" 25 cm



MUSEUX UTERINE
TENACULUM 24 CM



JACOBS UTERINE
TENACULUM 21 CM



GY-008-27
PRATT 10" 25 cm

UTERINE HOLDING FORCEPS/ GYNAECOLOGY



GY-009-27 COLLIN Uterine elevating Forceps 10", 25 cm




OB-003-21 Green-Armytage Caesarian Forceps 8 1/4"



UTERINE MANIPULATOR "HULKA"



SOMER Uterine seizing Forceps 22.5 cm 9"



UG-002-27 DARTIGUES Uterine seizing Forceps 28 cm 11"

Uro Gynaecology Retractors - Pelvic Surgery



RT-050-01 ST. MARKS
Pelvic Retractor



RT-058-28
SIMON Retractor



RT-059-28
BURNER Retractor



GY-024-28 LANDAU
Vaginal Specula
100X28 MM



GY-025-25 BREISKY
Vaginal
Specula 100X15 MM

Pelvic Surgery Ligature Carriers & Fascialata Stripper



“numed” MIYA ZAKI LIGATURE CARRIER FOR SACROSPINOUS SUSPENSION

Miya Hook Designed specifically for Prolapsed Vaginal Vault Repair

prolapsed vaginal vault repair the Miya Hook Ligature Carrier Set facilitates maximum exposure thus permitting a secure, reliable ligature placement.


- Pulls Down into the safety of Pararectal Space
- Helps avoid complication with prudenal
- Sciatic nerves and inferior gluteal vessels
- Requires less dissection and minimizing blood loss.



Suture Notch allows easier suture retrieval



Suture wedge anchors and organizes sutures



Nichols-Veronikis Ligature Carrier™ jaw closed on the sacrospinous ligament. Nichols Hook retrieving suture

“numed” Ligature Carrier for Sacrospinous Colpopexy

CMP-007-23 HEANY Hystrectomy Clamp 1 Tooth Heavy Pattern Non-Slip Graduated



Retractor for Upper Abdomen & Bladder Self Retaining



RT-058-01 O'Sullivan O'Coonner Retractor Ring
with Upper Abdominal Blade 110x70 mm and Lower
Bladder Blade 110x70 mm mm Self Retaining



"numed" Double Bar Self Retaining Abdominal
Retractor with Lower Bladder Blade and Upper
Abdominal Blade extended through 'C' Bar
Attachment

Introducing "nu" clamps

HYSTRECTOMY CLAMPS GREATER HOLDING POWER

Minimal Tissue Trauma
Minimal Tissue Slippage
Lasting Product Durability
Highest Quality Standard
CNC Manufactured

Leading surgeons from all over the nation perceive "nu" clamps as the premier instrument line for gynecological procedures. "nu med" Hysterectomy Clamp brings the second generation of Hysterectomy Clamps and Scissors which is CNC manufactured for ultimate precision and safety.



Numed Adnexal Clamp

Delicate instruments especially suitable for transecting the adnexa. The Adnexal Clamps are a delicate version of the Parametrium Clamps.

16171 GY 190 mm slightly curved
16172 GY 190 mm strongly curved
16271 GY 210 mm slightly curved
16272 GY 210 mm strongly curved



Numed Lateral Clamps

Clamps specially developed for resection of the lateral parametrium with radical hysterectomy (Wertheim). The S-shaped curvature is very easy to use and space saving

16577 GY 240 mm S-shaped
16578 GY 300 mm S-shaped

HYSTRECTOMY CLAMPS



Designed in response to modern biomechanical findings, these proven forceps afford a secure, non-slip hold with minimal tissue trauma. With its different overall and jaw lengths, as well as anatomically bent jaws. The Parametrium Clamps meet the high demands of a safe, non-slip and non crushing hysterectomy clamp.

Numed Parametrium Clamps

- 16070 GY** 210 mm straight
- 16071 GY** 210 mm slightly curved
- 16072 GY** 210 mm strongly curved
- 16073 GY** 210 mm angled



A universal instrument for abdominal and vaginal hysterectomy. The S-shaped shank and jaw design allows excellent vision into the Narrow operating field.

Numed Fischer Clamp

- 16471 GY** 220 mm slightly curved

HYSTRECTOMY CLAMPS



- 16175 GY** 240 mm straight
- 16176 GY** 240 mm slightly curved
- 16177 GY** 240 mm strongly curved
- 16178 GY** 240 mm angled

Numed Werthim-Cullin Vaginal Clamps

These clamps are especially suitable for vaginal Surgery. The proven jaw profile provides a secure, non-slip hold with minimal tissue trauma.



Numed Lainz Clamps

Infundibular Pelvic Clamps

The IFP (Infundibular Pelvic Clamp) is perfectly shaped to reach high into the pelvis during vaginal hysterectomy salpingo-oophorectomy (BSO). The s-shaped shanks of the instrument out of the way of the surgeon.

The 70mm jaw length is able to envelop the infundibular pelvic ligament and round ligament simultaneously The IFP-XL clamp offers the extra length of two inches proximal to the box lock to accommodate the larger patient.

- 16079 GY** 240 mm
- 16379 GY** 300 mm

Cardinal Ligament Clamp



16122 GY 290 mm Numed Cardinal Clamp

Especially “soft” cardinal clamp for resection of the highly vascular and delicate cardinal ligament



GY-13-0-11 Numed Multi Ligation Clamp 260 mm

For blunt vessel dissection in gynaecological vessels or vascular structures during ligation. This clamp is supplied with DeBakey jaw profile

Deaver Retractors For Abdominal & Vaginal Surgery

B1: "S" type Deaver retractors



B1-01: blade width 50mm
blade depth 155mm

B1-02: blade width 40 mm
blade depth 130mm

B1-03: blade width 22 mm
blade depth 115mm

B1-04: width 20mm
depth 210mm

B1-05: width 28mm
depth 230mm

B1-06: width 30mm
depth 280mm



Suction With Light

D4 Suction tube with light curved

D4-01 2.5

D4-02 3.0

D4-03 3.5

“L” Shape Retractors

B2: Retractor (Both sides with lip)



B3: “L” type Retractor



- B2-01: blade width 60 mm
blade depth 150mm
- B2-02: blade width 50 mm
blade depth 110mm
- B2-03: blade width 40 mm
blade depth 85mm

- B3-01: blade width 60 mm
blade depth 168mm
- B3-02: blade width 50 mm
blade depth 150mm
- B3-03: blade width 38 mm
blade depth 136mm
- B3-04: blade width 24 mm
blade depth 100mm

“L” Shape Retractors



B4: Retractor with light
Vertebra retractor

B4-01: width * depth
32mm*105mm

B4-02: 30mm*85mm

B4-03: 28mm*65mm

B5: Retractor with light

B5-01: width * depth
26mm*140mm

B5-02: 26mm*100mm

B5-03: 26mm*80mm

Anal Retractor



C2: Vaginal Retractor with light

(with two blades)

C2-01 width 22mm

depth 85mm

C2-02 width 22mm

depth 80mm

C2-03 width 22mm

depth 70mm



E3 Bladder retractors with light
Narrow Blades (With two
blades)

E3-01 width 35mm
depth 100mm

E3-02 : width 35mm
depth 105mm

L Shape Vaginal Speculum with light



C3-01: Vaginal retractor with light

width 20mm

depth 85mm

C3-02: Vaginal retractor with light

width 20mm

depth 85mm



E4: Gynecology retractor (arc type)

Pelvic Retractor with Light

Width: 21mm

Depth: 130mm